

# Animus<sup>®</sup> Monovite<sup>®</sup>

**Poly(glycolide cocaprolactone)  
(PGCL) – monofilament  
absorbable sterile surgical suture.**

## INDICATIONS

PGCL sutures are indicated for use in soft tissue approximation and/or ligation.

## CONTRAINDICATIONS

Not for use in cardiovascular, neurological surgery, microsurgery or ophthalmic surgery.

- Monovite has an excellent tensile strength with a reliable absorption rate of 91 – 120 days.
- Monovite has superb handling with exceptional suppleness and good knot security.

<b>70%</b>	tensile strength after	<b>7</b> DAYS
<b>40%</b>	tensile strength after	<b>14</b> DAYS
<b>12%</b>	tensile strength after	<b>21</b> DAYS

- Monovite should be stored in its sealed foil pouch until ready for use.
- The product should then be used within 30 days and stored at room temperature in a dry environment.



**Animus<sup>®</sup>**  
**Monovite<sup>®</sup>**

**Animus Monovite (PGCL) Cassettes**

ANIMUS CODE	DESCRIPTION	USP	METRIC	LENGTH	PRICE
AT2025	Animus Monovite Cassette	3-0	2	25m	45.75
AT3025	Animus Monovite Cassette	2-0	3	25m	47.75
AT3525	Animus Monovite Cassette	0	3.5	25m	56.25
AT4025	Animus Monovite Cassette	1	4	25m	62.50

**Animus Monovite (PGCL) Swaged (Packs of 10 – racetrack format)**

ANIMUS CODE	DESCRIPTION	USP	METRIC	LENGTH	PRICE
AT8201	Animus Monovite Reverse Cutting 24mm 3/8 circle	3-0	2	75cm	22.65
AT8202	Animus Monovite Taper Point 26mm 1/2 circle	3-0	2	75cm	23.35
AT8301	Animus Monovite Reverse Cutting 24mm 3/8 circle	2-0	3	75cm	24.35
AT8302	Animus Monovite Taper Point 26mm 1/2 circle	2-0	3	75cm	23.35