# ANIMUS SURGICAL LTD PRODUCT BROCHURE PRICE LIST & ORDER FORM

2023

P O Box 1260 Yateley GU47 7GH

T: 01252 845988

E: info@animussurgical.com





**Veterinary Suture Specialists** 

Violine Monovite Catgut Nylomide NyloBlue

# Animus® Catgut

Veterinary Suture Specialists



# Chromic catgut – absorbable sterile surgical suture composed of connective tissue.

#### **INDICATIONS**

Chromic gut suture is indicated for the use in animals in general soft tissue approximation and/or ligation, including use in ophthalmic surgery.

#### **CONTRAINDICATIONS**

Not for use in cardiovascular or neurological surgery.

- Purity of collagen for optimal tensile strength and in vivo performance.
- Surface smoothness to facilitate knot run down and knot strength.
- Absorption rate is variable at 14 28 days.
- Catgut should be stored sealed until ready for use.
- The product should be stored at room temperature in a dry environment.

### Animus Catgut

#### **Animus Catgut Cassettes**

ANIMUS CODE	DESCRIPTION	USP	METRIC	LENGTH	PRICE
AK260	Animus Catgut Cassette	4-0	2	50m	80.97
AK261	Animus Catgut Cassette	3-0	3	50m	68.11
AK262	Animus Catgut Cassette	2-0	3.5	50m	86.56
AK263	Animus Catgut Cassette	0	4	50m	94.58
AK264	Animus Catgut Cassette	1	5	50m	101.71
AK265	Animus Catgut Cassette	2	6	45m	112.80
AK266	Animus Catgut Cassette	3	7	40m	109.41
AK267	Animus Catgut Cassette	4	8	25m	81.20



**Veterinary Suture Specialists** 

Violine Monovite Catgut Nylomide NyloBlue

Animus®

### **Nylomide**

### NyloBlue

Nylomide – non-absorbable braided sterile surgical suture.


NyloBlue – non-absorbable monofilament sterile surgical suture.

#### **INDICATIONS**

Nylon sutures are indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological surgery.

#### CONTRAINDICATIONS

Due to the gradual loss of tensile strength that may occur over prolonged periods of in vivo, nylon sutures should not be used where permanent retention of tensile strength is required.



- Nylomide / NyloBlue should be stored sealed until ready for use.
- The product should be stored at room temperature in a dry environment.



## Animus<sup>®</sup> NyloBlue

#### **Animus Nylomide Cassettes**

ANIMUS CODE	DESCRIPTION	USP	METRIC	LENGTH	PRICE
AK271	Animus Nylomide Cassette	3-0	2	100m	38.44
AK272	Animus Nylomide Cassette	2-0	3	100m	40.18
AK275	Animus Nylomide Cassette	1	4	75m	40.18
AK273	Animus Nylomide Cassette	2	5	50m	49.50
AK276	Animus Nylomide Cassette	3+4	6	50m	49.50

#### **Animus NyloBlue Cassettes**

ANIMUS CODE	DESCRIPTION	USP	METRIC	LENGTH	PRICE
AB271	Animus NyloBlue Cassette	3-0	2	100m	40.92
AB272	Animus NyloBlue Cassette	2-0	3	100m	40.92



# Animus<sup>®</sup> Monoline<sup>®</sup>


Polydioxanone (PDO) – monofilament absorbable sterile surgical suture.



PDO sutures are intended for use in animals in soft tissue approximation and/or ligation, including use in paediatric cardiovascular and ophthalmic surgery.



Not for use in adult cardiovascular, microsurgery or neurological surgery.



• Monoline has superb handling with exceptional suppleness and excellent knot security.

75% tensile strength at
66% tensile strength at
53% tensile strength at
42 DAYS

- Monoline should be stored in its sealed foil pouch until ready for use.
- The product should then be used within 60 days and stored at room temperature in a dry environment.





### Animus<sup>®</sup> Monoline<sup>®</sup>

#### **Animus Monoline (PDO) Cassettes**

ANIMUS CODE	DESCRIPTION	USP	METRIC	LENGTH	PRICE
AM2015	Animus Monoline Cassette	3-0	2	15m	31.82
AM2030	Animus Monoline Cassette	3-0	2	30m	60.11
AM3015	Animus Monoline Cassette	2-0	3	15m	32.56
AM3030	Animus Monoline Cassette	2-0	3	30m	60.11
AM3515	Animus Monoline Cassette	0	3.5	15m	37.48
AM3530	Animus Monoline Cassette	0	3.5	30m	66.22
AM4015	Animus Monoline Cassette	1	4	15m	41.40
AM4030	Animus Monoline Cassette	1	4	30m	73.93

#### Animus Monoline (PDO) Swaged (Packs of 10 – racetrack format)

ANIMUS CODE	DESCRIPTION	USP	METRIC	LENGTH	PRICE
AM8151	Animus Monoline Reverse Cutting 19mm 3/8 circle	4-0	1.5	75cm	42.84
AM8152	Animus Monoline Taper Point 22mm 1/2 circle	4-0	1.5	75cm	42.84
AM8201	Animus Monoline Reverse Cutting 24mm 3/8 circle	3-0	2	75cm	45.43
AM8202	Animus Monoline Taper Point 26mm 1/2 circle	3-0	2	75cm	39.39
AM8301	Animus Monoline Reverse Cutting 24mm 3/8 circle	2-0	3	75cm	45.43
AM8302	Animus Monoline Taper Point 26mm 1/2 circle	2-0	3	75cm	45.43
AM8351	Animus Monoline Reverse Cutting 36mm 1/2 circle	0	3.5	75cm	45.43
AM8401	Animus Monoline Reverse Cutting 36mm 1/2 circle	1	4	75cm	39.39



# Animus® Violine®


# Polyglycolic acid (PGA) – braided absorbable sterile surgical suture.

#### INDICATIONS

PGA sutures are intended for use in animals in soft tissue approximation and/or ligation, including use in ophthalmic surgery.

#### **CONTRAINDICATIONS**

Not for use in cardiovascular or neurological surgery.



- Violine has superb handling with exceptional suppleness and knot security.
- Violine + has the added assurance of antibacterial protection.

83% tensile strength after
52% tensile strength after
21 DAYS
20% tensile strength after
28 DAYS

- Violine should be stored in its sealed foil pouch until ready for use.
- The product should then be used within 90 days and stored at room temperature in a dry environment.





#### **Animus Violine (PGA) Cassettes**

ANIMUS CODE	DESCRIPTION	USP	METRIC	LENGTH	PRICE
AV1515	Animus Violine Cassette	4-0	1.5	15m	29.60
AV1530	Animus Violine Cassette	4-0	1.5	30m	50.78
AV2015	Animus Violine Cassette	3-0	2	15m	29.60
AV2030	Animus Violine Cassette	3-0	2	30m	50.78
AV3015	Animus Violine Cassette	2-0	3	15m	33.06
AV3030	Animus Violine Cassette	2-0	3	30m	50.78
AV3515	Animus Violine Cassette	0	3.5	15m	33.06
AV3530	Animus Violine Cassette	0	3.5	30m	50.78
AV4015	Animus Violine Cassette	1	4	15m	33.06
AV4030	Animus Violine Cassette	1	4	30m	50.78
AV5015	Animus Violine Cassette	2	5	15m	33.93
AV6015	Animus Violine Cassette	3+4	6	15m	36.44

#### **Animus Violine + (PGA) Cassettes**

ANIMUS CODE	DESCRIPTION	USP	METRIC	LENGTH	PRICE
AP2025	Animus Violine + Cassette	3-0	2	25m	46.61
AP3025	Animus Violine + Cassette	2-0	3	25m	59.47
AP3525	Animus Violine + Cassette	0	3.5	25m	63.32
AP4025	Animus Violine + Cassette	1	4	25m	69.75

#### Animus Violine (PGA) Swaged (Packs of 10 – racetrack format)

ANIMUS CODE	DESCRIPTION	USP	METRIC	LENGTH	PRICE
AV8151	Animus Violine Reverse Cutting 19mm 3/8 circle	4-0	1.5	75cm	33.14
AV8152	Animus Violine Taper Point 22mm 1/2 circle	4-0	1.5	75cm	30.69
AV8201	Animus Violine Reverse Cutting 24mm 3/8 circle	3-0	2	75cm	33.14
AV8202	Animus Violine Taper Point 26mm 1/2 circle	3-0	2	75cm	30.69
AV8301	Animus Violine Reverse Cutting 24mm 3/8 circle	2-0	3	75cm	33.14
AV8302	Animus Violine Taper Point 26mm 1/2 circle	2-0	3	75cm	30.69
AV8351	Animus Violine Reverse Cutting 36mm 1/2 circle	0	3.5	75cm	33.14
AV8401	Animus Violine Reverse Cutting 36mm 1/2 circle	1	4	75cm	33.14
AV8402	Animus Violine Taper Point 37mm 1/2 circle	1	4	75cm	30.69



**Animus**<sup>®</sup>

**Monovite** 

# Animus<sup>®</sup> Monovite<sup>®</sup>

Poly(glycolide cocaprolactone) (PGCL) – monofilament absorbable sterile surgical suture.



PGCL sutures are indicated for use in soft tissue approximation and/or ligation.

#### CONTRAINDICATIONS

Not for use in cardiovascular, neurological surgery, microsurgery or ophthalmic surgery.



METRIC:

- Monovite has an excellent tensile strength with a reliable absorption rate of 91 120 days.
- Monovite has superb handling with exceptional suppleness and good knot security.

70% tensile strength after
40% tensile strength after
14 DAYS
12% tensile strength after
21 DAYS

- Monovite should be stored in its sealed foil pouch until ready for use.
- The product should then be used within 30 days and stored at room temperature in a dry environment.

### Animus<sup>®</sup> Monovite<sup>®</sup>

#### **Animus Monovite (PGCL) Cassettes**

ANIMUS CODE	DESCRIPTION	USP	METRIC	LENGTH	PRICE
AT2025	Animus Monovite Cassette	3-0	2	25m	61.56
AT3025	Animus Monovite Cassette	2-0	3	25m	64.25
AT3525	Animus Monovite Cassette	0	3.5	25m	74.00

#### Animus Monovite (PGCL) Swaged (Packs of 10 – racetrack format)

ANIMUS CODE	DESCRIPTION	USP	METRIC	LENGTH	PRICE
AT8201	Animus Monovite Reverse Cutting 24mm 3/8 circle	3-0	2	75cm	37.09
AT8202	Animus Monovite Taper Point 26mm 1/2 circle	3-0	2	75cm	39.68
AT8301	Animus Monovite Reverse Cutting 24mm 3/8 circle	2-0	3	75cm	39.68
AT8302	Animus Monovite Taper Point 26mm 1/2 circle	2-0	3	75cm	32.20





**Membra** is a long-lasting, transparent bandage replacement made of a specially formulated, non-toxic, biodegradable hydrogel which dries to form a water-tight dressing.

It is designed for the protection of open wounds and to strengthen surgical incisions.

Membra can be used across a range of species (including aquatic animals) and one application lasts up to 14 days.

Hydrogels, such as PEG (Polyethylene glycol) offer advantages compared to other wound dressings:

- hydration of the wound bed promote healing and autolytic debridement of necrotic tissue.
- cooling of the wound surface –soothe and are a good choice to protect skin affected by burning.
- transparent healing can be monitored without the need for dressing changes.
- non-toxic and non-irritant won't encourage licking or chewing of affected areas.

When applied to a dry wound, Membra dries in a few minutes to form a water- resistant seal.

- Ideal for animals that live outside, in water or are exposed to the elements.
- In most cases there is no need for bandaging saving time, cost and owner inconvenience.

#### Membra promotes healing:

- Allows oxygen through to the wound surface.
- · Blocks external contaminants.
- Maintains the ideal moisture levels for healing.

Membra is durable, flexible and biodegrades naturally.

- Lasts up to 14 days.
- Can be used on hard-to-bandage areas.
- Allows the animal to move and return to normal function.
- Doesn't require removal.

#### In a wide range of species

- > All small mammals.
- Equine and livestock species.
- Avian species.
- Exotics and reptiles e.g. snakes, lizards, frogs.
- Aquatic and marine species including fish, sea turtles, whales and seals.
- Zoo Animals.

#### For a wide range of lesions

- ✓ For added protection on top of sutures lines.
  - Orthopaedic procedures.
  - Reconstructive surgery.
  - Lump / tumour removal.
  - Amputations.
- ✓ To reduce / eliminate bandaging for open wounds.
  - Trauma/Emergency wounds.
  - Degloving injuries.
  - Following large mass removal.
- ✓ Lacerations / scrapes / cuts / abrasions.
- ✓ Burns.
- ✓ Lick granulomas.
- ✓ Paw pad cuts and tears.
- ✓ Ideal for hard-to-bandage locations.

Animus Code: AS0102 - Price - £69 excl VAT

### ANIMUS MEDBOND SKIN GLUE Cyanoacrylate Skin Adhesive A quick, painless solution that helps prevent infection



#### **DESCRIPTION:**

Cyanoacrylate Skin Adhesive

#### **FEATURES:**

- Enables precise placement of adhesive.
- · Adhesive stops minor bleeding.
- Acts as a barrier to help keep foreign matter from entering the wound.
- May eliminate need for bandages in some cases.
- Saves time and money for minor procedures.

#### **INDICATIONS**

Animus Medbond Skin Glue liquid adhesive bandage is intended for surgical incisions, lacerations, pad wounds and bleeding nails. Never apply ANY type of adhesive IN a wound. Intended application is at the edges of the area to be treated. Do not apply Animus Medbond Skin Glue to the eyes. Do not use on infected areas, or wounds that are draining. Never use on animal bites. Do not use on mucosal surfaces (eg, oral cavity, lips).

#### COLOR(S)

Violet

#### **PACKAGING**

Included in each pack: 1 Animus Medbond Skin Glue bottle and 12 non-clogging applicators.

FOR VETERINARY USE ONLY

Animus Code: AG0102 Price £13.02 excl VAT

T 01252 845988 E info@animussurgical.com www.animussurgical.com

ANIMUS SURGICA	AL LTD PRODUCTS LIST & ORDER FORM					
ANIMUS CATGUT		USP	Metric	Length	RRP	No:
	Animus Catgut Cassette	4-0	2	50m	80.97	NO.
	Animus Catgut Cassette	3-0	3	50m	68.11	
	Animus Catgut Cassette	2-0	3.5	50m	86.56	
	Animus Catgut Cassette	0	4	50m	94.58	
	Animus Catgut Cassette	1	5	50m	101.71	
	Animus Catgut Cassette	2	6	45m	112.80	
	Animus Catgut Cassette	3	7	40m	109.41	
	Animus Catgut Cassette	4	8	25m	81.20	
ANIMUS NYLON C		-			01120	
	Animus Nylomide Cassette	3-0	2	100m	38.44	
	Animus Nylomide Cassette	2-0	3	100m	40.18	
	Animus Nylomide Cassette	1	4	75m	40.18	
	Animus Nylomide Cassette	2	5	50m	49.50	
	Animus Nylomide Cassette	3+4	6	50m	49.50	
	Animus NyloBlue Cassette	3-0	2	100m	40.92	
	Animus NyloBlue Cassette	2-0	3	100m	40.92	
	(PGA) CASSETTES			100111	40102	
	Animus Violine Cassette	4-0	1.5	15m	29.60	
	Animus Violine Cassette	4-0	1.5	30m	50.78	
	Animus Violine Cassette	3-0	2	15m	29.60	
	Animus Violine Cassette	3-0	2	30m	50.78	
	Animus Violine Cassette	2-0	3	15m	33.06	
	Animus Violine Cassette	2-0	3	30m	50.78	
	Animus Violine Cassette	0	3.5	15m	33.06	
	Animus Violine Cassette	0	3.5	30m	50.78	
	Animus Violine Cassette	1	4	15m	33.06	
	Animus Violine Cassette	1	4	30m	50.78	
	Animus Violine Cassette	2	5	15m	33.93	
	Animus Violine Cassette	3+4	6	15m	36.44	
	+ (PGA) CASSETTES					
	Animus Violine + Cassette	3-0	2	25m	46.61	
	Animus Violine + Cassette	2-0	3	25m	59.47	
	Animus Violine + Cassette	0	3.5	25m	63.32	
	Animus Violine + Cassette	1	4	25m	69.75	
	NE (PDO) CASSETTES			-		
	Animus Monoline Cassette	3-0	2	15m	31.82	
	Animus Monoline Cassette	3-0	2	30m	60.11	
	Animus Monoline Cassette	2-0	3	15m	32.56	
	Animus Monoline Cassette	2-0	3	30m	60.11	
	Animus Monoline Cassette	0	3.5	15m	37.48	
	Animus Monoline Cassette	0	3.5	30m	66.22	
	Animus Monoline Cassette	1	4	15m	41.40	
	Animus Monoline Cassette	1	4	30m	73.93	
	TE (PGCL) CASSETTE					
	Animus Monovite Cassette	3-0	2	25m	61.56	
	Animus Monovite Cassette	2-0	3	25m	64.25	
	Animus Monovite Cassette	0	3.5	25m	74.00	

		USP	Metric	Length	RRP	No:
<b>ANIMUS VIOLII</b>	NE (PGA) SWAGED PACKS OF 10 IN RACETRACK FORMAT					
AV8151	Animus Violine Reverse Cutting 19mm 3/8 circle	4-0	1.5	75cm	33.14	
AV8152	Animus Violine Taper Point 22mm 1/2 circle	4-0	1.5	75cm	30.69	
AV8201	Animus Violine Reverse Cutting 24mm 3/8 circle	3-0	2	75cm	33.14	
AV8202	Animus Violine Taper Point 26mm 1/2 circle	3-0	2	75cm	30.69	
AV8301	Animus Violine Reverse Cutting 24mm 3/8 circle	2-0	3	75cm	33.14	
AV8302	Animus Violine Taper Point 26mm 1/2 circle	2-0	3	75cm	30.69	
AV8351	Animus Violine Reverse Cutting 36mm 1/2 circle	0	3.5	75cm	33.14	
AV8401	Animus Violine Reverse Cutting 36mm 1/2 circle	1	4	75cm	33.14	
AV8402	Animus Violine Taper Point 37mm 1/2 circle	1	4	75cm	30.69	
ANIMUS MONO	LINE (PDO) SWAGED PACKS OF 10 IN RACETRACK FORMAT					
AM8151	Animus Monoline Reverse Cutting 19mm 3/8 circle	4-0	1.5	75cm	42.84	
AM8152	Animus Monoline Taper Point 22mm 1/2 circle	4-0	1.5	75cm	42.84	
AM8201	Animus Monoline Reverse Cutting 24mm 3/8 circle	3-0	2	75cm	45.43	
AM8202	Animus Monoline Taper Point 26mm 1/2 circle	3-0	2	75cm	39.39	
AM8301	Animus Monoline Reverse Cutting 24mm 3/8 circle	2-0	3	75cm	45.43	
AM8302	Animus Monoline Taper Point 26mm 1/2 circle	2-0	3	75cm	45.43	
AM8351	Animus Monoline Reverse Cutting 36mm 1/2 circle	0	3.5	75cm	45.43	
AM8401	Animus Monoline Reverse Cutting 36mm 1/2 circle	1	4	75cm	39.39	
ANIMUS MONO	VITE (PGCL) SWAGED [PACKS OF 10 IN RACETRACK FORMAT					
AT8201	Animus Monovite Reverse Cutting 24mm 3/8 circle	3-0	2	75cm	37.09	
AT8202	Animus Monovite Taper Point 26mm 1/2 circle	3-0	2	75cm	39.68	
AT8301	Animus Monovite Reverse Cutting 24mm 3/8 circle	2-0	3	75cm	39.68	
AT8302	Animus Monovite Taper Point 26mm 1/2 circle	2-0	3	75cm	32.20	
ANIMUS MEDB	OND SKIN GLUE					
AG0102	Animus Medbond Skin Glue				13.02	
ANIMUS MEME	RA WOUND SEALENT					
AS0102	Animus Membra Wound Sealant				69.00	

Practice Name & Addi	ıs – <u>info@anımussu</u>	rgical.com
		••••••
Email address:	 	



Violine Monoline Monovite Catgut Nylomide NyloBlue